

План статьи

Введение

1. Какую информацию об образце мы можем получить, используя метод ЯМР в низких полях?
2. О молекулярных движениях, как источнике информации, с точки зрения ЯМР-релаксометрии. О видах движения (работы Кая).
3. О влиянии времени корреляции на T_1 , T_2 твердых и жидких фаз образцов и форму спектральной линии.
4. Корреляция структуры и свойств веществ с амплитудными и временными характеристиками спадов и параметров спектра.
5. Микрообзор того, что кто смотрит по фиду.
6. ГЭП – разброд и шатание в подходах к обработке фидов.
7. Цель работы.

Результаты и Дискуссия

1. Твердые тела. Подбор ЯМР-экспериментов с кратким описанием их особенностей.
2. Обработка данных во временной области.
3. Модели Абрагам, Гаусс, Вейбулл, Эксп. в зависимости от молекулярной динамики. Комбинации моделей в гетерофазных образцах. Модели не всегда идеально подходят – наши экспериментальные данные по разным видам твердых образцов.
4. Обработка данных в частотной области.
5. Формы спектральных линий. Второй момент. Взаимосвязь второго момента с временем поперечной релаксации.
6. Новый метод: Second moment approximation.
7. Идея.
8. Примеры применения к спадам с осцилляцией и без.
9. Сопоставление с результатами подходов Интегрирования и Усреднения.

10. Возможные будущие задачи. Может, как-то корреляцию с DQ поизучать.
11. Медленные движения молекул твердых тел эффективнее исследовать по спаду свободной индукции.
12. Есть модификации как RK, Solid Echo (самые сильные дипольные взаимодействия), MSE (вроде, слабые).
13. На каждом из упомянутых экспериментов по-разному отражается константа остаточных дипольных взаимодействий.

Analysis of FID shapes of Solids in glassy condition

Baranauskaite Valeriia, Mariia Ivanova, Leonid's Sister, Leonid Grunin

23 октября 2024 г.

Abstract

Why is it important

What has been done

What was found

What is the general result

How it can be used in the future

1 Introduction

2 Experiment details

3 Results and discussion

3.1 FID processing in Time-Domain NMR

The recorded FID signal in NMR experiments typically consists of time, and real and imaginary components. Before analyzing the shape of FID, it undergoes preprocessing steps, including baseline subtraction, phase and frequency adjustments in the time domain. The NMR signal for baseline is recorded in the empty ampule at exactly the same conditions as the sample. Phase adjustment addresses phase errors in the time domain data. These errors can arise from imperfections in the experimental setup or sample characteristics. The process involves finding the optimal phase angle that minimizes the difference between the real component of the FID signal and the overall amplitude of the signal within the initial part (typically before 30 microseconds). By adjusting the phase angle, we aim to align the real component of the signal with its amplitude, ensuring accurate representation of the signal's characteristics.

Frequency adjustment ensures that the frequency spectrum of the signal is correctly positioned and scaled. This adjustment is crucial due to potential imperfections in the equipment or sample, which can cause the signal to deviate from its ideal state. The Fast Fourier Transform (FFT) procedure transforms the normalized time-domain signal into a frequency-domain in the form of a spectrum. The frequency axis is determined based on the time axis, following the Nyquist theorem. The spectrum is then shifted so that the maximum aligns with zero frequency, ensuring an even representation of positive and negative frequencies, i.e. symmetric spectrum. Finally, an inverse FFT (iFFT) is performed to obtain the FID signal with the corrected frequency.




Рис. 3.1: The comparison between original (solid lines) and corrected (dashed lines) signals

To mitigate the effects of magnet inhomogeneity at short time scales (before $30\mu s$), where

the fast relaxation occurs, it is recommended to subtract the long component. Glycerol, known for its long relaxation times, serves as a suitable reference. Alternatively, any substance with a transverse relaxation time exceeding 2 milliseconds can be used. The Glycerol's reference FID is normalized, adjusted for frequency and phase in the same manner as it was done with Sample. Subsequently, the real component of the Sample is divided by the amplitude of the Glycerol FID.

To eliminate the long component, an exponential fitting of the real part of the FID corresponding to the slow decay within the time range before $200\mu s$ should be conducted and the fitted curve should be subtracted from the original data.




Рис. 3.2: The comparison of the original FID (red) and subtracted FID (blue)

The concluding steps in the time-domain involve apodizing the real and imaginary components of the FID, followed by zero-filling to ensure that the final FID comprises a power of 2 number of data points (2^n). This step is crucial as FFT computation is more efficient with data containing this specific number of points. Subsequently, the FFT procedure is applied, yielding the resulting spectrum in the frequency-domain.

$$f(\text{apodization}) = \exp\left(-\left[\frac{\text{Time}}{\sigma}\right]^4\right) \quad (1)$$

3.2 FID processing in Frequency-Domain NMR

The spectra derived from the FFT procedure can be manually phased using linear or polynomial functions. Subsequently, the resulting real and imaginary components of the spectra can be recalculated accordingly.

$$\begin{aligned}
Re &= Re \cdot \cos(\phi) - Im \cdot \sin(\phi) \\
Im &= Re \cdot \sin(\phi) + Im \cdot \cos(\phi) \\
Amplitude &= \sqrt{Re^2 + Im^2}
\end{aligned} \tag{2}$$

Additionally, apodization of the spectrum, akin to the time-domain signal, should be performed to ensure zero amplitude at high frequencies, utilizing similar equations as in Eq. 1. Once apodized, the acquired spectrum is prepared for further analysis, such as the calculation of the second moment.

4 SE and FID

In the analysis of Solid Echo (SE) and Free Induction Decay (FID) signals, several steps are undertaken to ensure precise data acquisition, correction, and interpretation as is described in Avsnittene 3.1 og 3.2. The SE signals and FID are shown in Fig. 4.1 a). The maxima of the SE signals are then plotted versus corresponding echo times as shown in Fig. 4.1 b). To interpret the data, a Gaussian function is employed for fitting.

$$G(x) = A \cdot \exp\left(-\frac{x^2}{2\sigma^2}\right) \tag{3}$$

The Gaussian function is of the standard form, where the key parameter A represents the amplitude. This value is directly derived from the SE signals and offers an essential reference for later stages of the analysis.




Рис. 4.1: a: The overall representation of SE signals and FID (red); b: The SE maxima versus the corresponding echo times

Next, the initial portion of the FID signal, particularly the segment between 10 and 20 microseconds is fitted using the same Gaussian model employed earlier for the SE Eq. 3. It is crucial to ensure that the amplitude of the Gaussian fit matches the maximum amplitude observed at time $t = 0$, as derived from the SE signal maxima. After fitting the FID signal, intersection points between the fitted Gaussian curve and the original FID signal are determined. These intersection points highlight critical aspects of the experimental data and assist in reconstructing the hidden portions of the FID using the fitted Gaussian model Fig. 4.2.




Рис. 4.2: The rebuild FID and the original FID

The final step involves validating the amplitude of the FID, which is rebuilt from the SE maxima. To ensure the accuracy of this reconstruction, the FID of a water sample with a known mass is recorded under identical experimental conditions. The maximum amplitudes of the FIDs for both the water and the cellulose samples are then compared, as these amplitudes correspond directly to the number of protons in each sample. The proton density, calculated using the formula:

$$\text{Proton density} = \frac{\max(\text{Amplitude})}{\left(\frac{\text{mass}}{M}\right) \cdot N_A \cdot (\text{number of protons per molecule})} \quad (4)$$

This formula provides a reference value for the proton density. By multiplying the proton density obtained from the water sample by the number of protons in the cellulose sample, a reference value of 42 is obtained. This reference is then compared to the maximum amplitude of 44 from the rebuilt FID. The close agreement between these values, within the margin of error, confirms the reliability of the proposed method for reconstructing the FID from SE signals.




Рис. 4.3: The recorded FID of water and cellulose sample

Compound	Mass (g)	Molar Mass (g/mol)	Proton Density (e-21)	FID Max Amplitude
Water	0.0963	18.01528	8.522	55
Cellulose	0.1334	162.1406	8.949	44

Таблица 1: Physical properties and FID amplitudes of Water and Cellulose samples

5 Conclusions

6 References